

Evaluation of the Formaldehyde Hemiacetals and Acetals Relevant to Electronic Cigarettes

Xiaohong Cathy Jin, Karen C. Avery, Regina M. Ballentine, William P. Gardner, Willie J. McKinney, Matt S. Melvin, Yezdi B. Pithawalla, Donna C. Smith and Karl A. Wagner



Altria

Altria Client Services


Background

- Formaldehyde is recommended for reporting in the FDA draft guidance for e-vapor products¹
- Methods for determination of formaldehyde typically use 2,4-Dinitrophenylhydrazine (DNPH) derivatization followed by analysis of the resulting hydrazone:
 - EPA method 8315A for waste water, air samples, etc.²
 - CORESTA CRM 74 for cigarette smoke³
 - Various methods for e-vapor products⁴



“Hidden Formaldehyde” (Jensen *et al.*)

- Reported presence of formaldehyde hemiacetals of propylene glycol (PG) or glycerin (Gly) in e-vapor aerosol⁵
- Suggested conversion of hemiacetal to acetal in acidic DNPH solution⁶



- Claimed that the commonly-used DNPH methods underestimate formaldehyde levels in e-vapor




Objective

- Evaluate the DNPH method performance for analysis of formaldehyde in e-vapor aerosol:
 1. Do PG- and Gly-hemiacetals hydrolyze to formaldehyde in acidic DNPH solution?
 2. Do PG- and Gly-hemiacetals convert to the corresponding acetals in acidic DNPH solution?
- Evaluate if acetals are formed in e-vapor aerosol




Do PG and Gly-hemiacetals Hydrolyze to Formaldehyde in Acidic DNPH Solution?




Approach

- Do hemiacetals hydrolyze to formaldehyde?



PG-hemiacetal (PG-HA)*



Gly-hemiacetal (Gly-HA)*

- Add PG-HA or Gly-HA to acidic DNPH solution
- Determine formaldehyde by UPLC/MS

* Custom synthesized by Dr. Sönke Peterson (Worms, Germany)



LC/MS Method for the Analysis of Formaldehyde

- Do hemiacetals hydrolyze to formaldehyde?

Parameters	
Instrumentation	UHPLC-MS
Column	BEH C18 (1.7 μ m)
Ionization mode	ESI negative, SIM
Trapping solution	DNPH in acetonitrile w/perchloric acid
Internal standard	Formaldehyde-DNPH-d3
Calibration range	0.01 - 4 μ g/mL (1.5 - 60 μ g/g)
Limit of quantification	1.5 μ g/g
Recovery	97 - 105%



Results

- Do hemiacetals hydrolyze to formaldehyde?

	Formaldehyde theoretical (µg/mL)	Formaldehyde determined (µg/mL)	Recovery (%)
PG-HA (3.19 µg/mL)	0.93	0.96	103
Gly-HA (5.86 µg/mL)	1.44	1.50	104


- PG-HA and Gly-HA quantitatively release formaldehyde and form the corresponding hydrazone in acidic DNPH solution
- These results corroborate those reported by PMI (Knorr, *et al.*)⁷




Summary

- Do hemiacetals hydrolyze to formaldehyde?

- PG-HA and Gly-HA quantitatively hydrolyze to formaldehyde in acidic DNPH solution
 - These hemiacetals rapidly hydrolyze in acidic DNPH solution to release formaldehyde
 - The secondary reaction with DNPH to form the hydrazone drives the equilibrium to release more formaldehyde




Do PG-HA and Gly-HA Convert to the Corresponding Acetals in Acidic DNPH Solution?




Approach

- Do hemiacetals convert to acetals?

- PG-HA and Gly-HA added to acidic DNPH solution ($\sim 5\mu\text{g/mL}$)
- Determine PG- and Gly-acetals by GC/MS



PG-acetal
(4-methyl-1,3-dioxolane)



Gly-acetals
(4-hydroxymethyl-1,3-dioxolane and
5-hydroxy-1,3-dioxane)



GC/MS Method for the Analysis of Acetal



- Do hemiacetals convert to acetals?

Parameters	
Instrument	GC/MS
GC column	Restek Rtx 624
Ionization mode	EI
Acquisition mode	SIM
Extraction / trapping	Dichloromethane (DCM)
Internal standard	2,3-Hexadione
Calibration range	0.01 - 1 $\mu\text{g/mL}$ (0.8 - 80 $\mu\text{g/g}$)
Limits of quantification	0.8 $\mu\text{g/g}$



Results

- Do hemiacetals convert to acetals?




➤ PG- and Gly-acetals were not formed



Summary

- Do hemiacetals convert to acetals?



- PG-HA and Gly-HA do not convert to the corresponding acetals in acidic DNPH solution



Are PG-Acetal or Gly-Acetals Formed in Aerosol?



Approach


- Do acetals form in e-vapor aerosol?

- Investigate acetal formation during aerosol generation by analyzing the e-liquid and the corresponding e-vapor



Experimental

- Do acetals form in e-vapor aerosol?



Results

- Do acetals form in e-vapor aerosol?

Sample	e-Liquid Composition	PG-Acetal (µg/g)		Gly-Acetal (µg/g)	
		e-liquid	aerosol	e-liquid	aerosol
Control	PG:Gly (50:50) (2.5% Nic, 15% H2O)	ND*	ND*	ND*	ND*
e-Cig A	Menthol (3.5% Nic)	ND*	ND*	<0.8	<0.8
e-Cig B	Non Menthol (3.5% Nic)	ND*	ND*	1.2±0.1	1.0±0.1

ND*: Not detected or <LOD (PG-acetal: 0.2ng/g, Gly-acetal: 0.3ng/g)

- PG- and Gly-acetals do not form during aerosolization process for the cig-a-like devices tested



Conclusions

- PG- and Gly-hemiacetals can be determined and reported as formaldehyde when using acidic DNPH
 - Quantitatively hydrolyze to release formaldehyde which forms the hydrazone in acidic DNPH solution
 - Do not convert to the respective acetals in acidic DNPH solution
- PG- and Gly-acetals are not formed during the aerosolization process for cig-a-like type products
 - Our results demonstrate that commonly-used DNPH methods are fit for the analysis of formaldehyde in e-vapor products
 - Our results **contradict** the reported claim that “Hidden Formaldehyde” in e-vapor products causes the underreporting of formaldehyde



References

1. FDA Draft Guidance for Industry, Premarket Tobacco Product Applications for Electronic Nicotine Delivery Systems, 2016
2. EPA method 8315A “Determination of Carbonyl Compounds by High Performance Liquid Chromatography” Rev.1, 1996
3. CORESTA CRM 74 “Determination of Selected Carbonyls in Mainstream Cigarette Smoke by HPLC”, Rev.4, 2018
4. Flora, J.W., Wilkinson, C.T.; Wilkson J. W.; Lipowicz P. J.; Skapars J. A.; Anderson, A.; and Miller J. H., Method for determination of Carbonyl Compounds in E-Cigarette Aerosols, J Chromatogr Sci. 2017 Feb; 55(2): 142-148
5. Jensen R.P.; Luo W.; Pankow J. F.; Strongin R. M.; and Peyton D. H., Hidden Formaldehyde in E-Cigarette Aerosols, N. ENGL. J. MED. 372:4, NEJM.ORG, JANUARY 22, 2015
6. Salamanca J.C.; Ian, M.; Escobedo J. O.; Jensen R. P.; Shaw A.; Campbell R.; Luo W.; Peyton D. H.; Strongin R. M., Formaldehyde Hemiacetal Sampling, Recovery, and Quantification from Electronic Cigarette Aerosols, Scientific reports, 7:11044, September 8, 2017,
7. Knorr A.; Gautier L.; Tekeste E.; Buchholz C., Almstertter M.; Arndt D.; Bently M., Formaldehyde-Glycerol Hemiacetal _Absence of “Hidden” Formaldehyde in THS 2.2 Aerosols, <https://www.researchgate.net/publication/291833243>



Reducing risk. Expanding choice.

Altria.

For copies of this presentation visit the Altria's Science
Website at www.altria.com/alcs-science

